



Green Lion

Digital Hearing Aid

SKU: GNINSAIDSBK
Model No: GL-HA2

Table of Contents

Overview	2
Specification	2
Application Scope	3
Touch Instructions	3
Hearing Aid Operation	3
Connecting to Wireless	4
Troubleshooting	5
Hearing Aid Adaptation Guide	6
Maintenance and Care	7
FCC Warning	7
Disposal	9
Warranty	9
Contact Us	9

Overview

The Green Lion Inspire Digital Hearing Aids provide advanced sound amplification with anti-whistle technology for enhanced hearing clarity. Offering up to 36 hours of battery life, these hearing aids come in a sleek charging case for easy storage and on-the-go convenience. Ideal for those looking to rediscover the joy of sound, the Inspire aids deliver a comfortable, discreet, and high-quality hearing experience.

Specifications

Model No.	GL-HA2
Material	ABS
Case Battery	650mAh
Earbuds Battery	50mAh
Max Sound Output	107 dB
Amplification Gain	Up to 30 dB
Noise Level (EIN)	≤25 dB
Volume Control	5 Stages
Bluetooth	V5.3
Transmission Range	10 m
Charging Port	USB-C
Charging Time	1.5 Hrs
Single Playtime	6 Hrs
Playtime With Case	36 Hrs
Product Weight	62 g
Product Dimensions	35.7*55.8*64 mm

Application Scope

This hearing aid is designed to aid patients with air-conduction hearing loss after professional fitting. It is not recommended for individuals with acute external otitis, tympanitis, otitis media, or those allergic to the materials used in the device.

The product uses a digital method for signal processing, converting the received sound signal (analog) into a digital signal. After processing and amplification, it is then converted back into an analog sound signal for output.

Touch Instructions

1. Play/Pause:

Long press the L/R.

2. Volume Adjustment:

2.1 Increase volume: Briefly press the L/R.

2.2 Decrease volume: Double-tap the L/R.

3. Calls

3.1 Answer a call: Briefly press the L/R.

3.2 Hang up a call: Long press the L/R.

3.3 Reject a call: Double-tap the L/R.

Note: These operations are only available when the device is connected via wireless.

Hearing Aid Operation

1. Take out the hearing aid from the charging case. The product will turn on automatically. After 15 seconds, the hearing aid will begin assisting hearing.

2. Volume and Mode Switching Operations table

Switch	Operation	Prompt tone
Volume level 1	Click the L/R	Minimum volume
Volume level 2	Click the L/R	level 2
Volume level 3	Click the L/R	level 3
Volume level 4	Click the L/R	level 4
Volume level 5	Click the L/R	Maximum volume
Strong noise cancellation mode	Double-tap the L/R	Strong noise cancellation mode
Indoor mood	Double-tap the L/R	Indoor mood
Outdoor mood	Double-tap the L/R	Outdoor mood

3. Turning off the device .To turn off the device, either place the earbuds back into the charging case, or long press the touch area of the earbuds outside the case for 6 seconds to turn them off individually.

Connecting to Wireless

The hearing aid can be connected to a mobile phone's Wireless feature for use as wireless earbuds without the need for hearing assistance. Follow the steps below:

- 1.** Take out the earphones and wait for Wireless pairing. The white/blue indicator light will flash alternately.
- 2.** Search for the Wireless name "GL-HA2" on the mobile phone's Wireless settings. After successful pairing, the white indicator light will flash slowly and then go out, with a prompt sound saying "connected."

Note:

- 1.** When the Wireless mode is connected and there is no sound output from the Wireless device, the hearing aid will automatically switch to the hearing aid mode within about 5 seconds.
- 2.** If sound is output from the Wireless device, the hearing aid will immediately exit hearing aid mode and switch to Wireless mode.

Troubleshooting

Issue	Cause & Resolution
Wireless connection issue (only one earbud connects or no connection)	<p>The earbuds may be connected to a previously paired device.</p> <p>Solution: Press and hold the RESET button on the charging case to clear previous connections, then re-pair the earbuds with the desired device.</p>
Squealing or noise	<p>Possible causes:</p> <ul style="list-style-type: none">- Interference from a strong magnetic field environment.- Improper earbud fit.- Incorrect earbud positioning. <p>Solution:</p> <ul style="list-style-type: none">- Move away from strong magnetic fields.- Ensure the earbud tips are the correct size and seal the ear canal properly.- Verify that the hearing aid is positioned correctly in the ear.
Weak or unclear sound	<p>Possible causes:</p> <ul style="list-style-type: none">- Low battery charge.- Improper earbud fit.- Foreign matter obstructing the earbud. <p>Solution:</p> <ul style="list-style-type: none">- Ensure the device is fully charged and try again.- Increase the volume.- Check for any debris blocking the earbud and ensure the earbud is securely placed in the ear.

Hearing Aid Adaptation Guide

Stage 1	Initial Familiarization Begin by acclimating to the sounds in your immediate environment, such as conversations, air conditioners, and running water from faucets. It is recommended to engage in conversations at a close distance to gradually adjust to the new hearing experience.
Stage 2	Outdoor Exploration Try walking in your neighborhood or a nearby park to observe how you perceive sounds from a distance. You may notice a clearer hearing experience for sounds that were previously less discernible. During this stage, it is recommended to wear your hearing aids for up to 4 hours daily for at least one week.
Stage 3	Engaging in Public Spaces You may begin to explore more populated environments such as local markets and shopping malls. Engage in activities where you can hear voices clearly, such as exercising outdoors or socializing. During this stage, aim to wear the hearing aids for at least 6 hours per day to avoid fatigue and ensure optimal adjustment.
Stage 4	Full Integration into Daily Life At this stage, you should be able to wear your hearing aids comfortably for 8 hours or more each day. You can fully immerse in complex sound environments like watching TV, listening to music, and going grocery shopping. If this is your first time wearing hearing aids, you may initially experience some discomfort or a sensation of fullness in the ears. Additionally, hearing your own amplified voice is a common adjustment. This is a natural process and may take 1 to 3 months for full adaptation.

Maintenance and Care

1. Cleaning the Earplugs


To clean the earplugs, carefully remove them from the hearing aids. Rinse with clean water or wipe them with alcohol wipes. After cleaning, make sure the earplugs are fully dried before reattaching them to the hearing aids.

2. Cleaning the Device

For the hearing aid device, use a lint-free soft cloth, dry cotton wool, or a soft sponge to wipe off any grease or moisture on the outer surface. Take extra care to prevent liquid from entering the device, including the charging cases.

3. Earplug Fitting

Users can select earplugs that best fit their individual ear shape (e.g., experience with ear canal size or if earplugs tend to fall out). Choose from the three available sizes: large (L), medium (M), and small (S). Once selected, attach the earplugs to the hearing aids at the sound outlet.



FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

1. Reorient or relocate the receiving antenna.
2. Increase the separation between the equipment and receiver.
3. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
4. Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirements. It can be used in portable exposure conditions without restriction.



Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

Products that you buy directly from our **Green Lion** website or shop come with a 24-month warranty.

When you buy **Green Lion** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at <https://www.greenlion.net/warranty> and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:
<https://www.greenlion.net/warranty>

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

Website: <https://www.greenlion.net/>

Sales: sales@greenlion.net

Service Support: Support@greenlion.net

Telegram: [@greenlion_global](https://t.me/greenlion_global)

Instagram: [@greenlion_global](https://www.instagram.com/greenlion_global)