



Green Lion

Black Head Remover

SKU: GNBKHD51RMWH

Model No: GL-TM51

Table of Contents

Preface	2
Safety Attention	2
Usage Prohibitions	3
Product Introduction	3
Features	4
Component Description	5
Beauty Suction Head	5
Replacement of Cosmetic Suction	7
Specifications	8
Disposal	9
Warranty	9
Contact Us	9

Preface

Thank you for purchasing the blackhead remover device. Please read the entire contents of this manual before use. This user manual contains important safety guidelines and instructions for the proper use of this product. For your safety and to prevent potential harm or damage from improper use, we strongly recommend that you carefully review the instructions. Keep this manual for future reference. This product is intended for personal use only and should not be used by multiple individuals in beauty salons or other professional settings.

Safety Attention

1. Please ensure that the filter sponge has been properly installed in the filter cup before use, and the sponge is dry, with no inhalation of water or other liquids.
2. At the same time, dry the water or other liquid on the skin where needed to be nursed, otherwise the dirty water or other liquid will be sucked into the air pump and cause blockage, preventing the machine from working properly.
3. The ornaments (such as necklaces, earrings, etc.) and contact lenses should be removed, otherwise, they could be sucked into the machine.
4. Please wipe the water or other liquid on your hands before nursing, to prevent the water or other liquid from entering the machine and causing damage to it.
5. This product can be used for the eyes, eyelids, ears, and neck, otherwise it may cause damage to the skin and body.
6. This product is not suitable for parts where eyes have taken cosmetic surgery, nasal augmentation, or other skin cosmetic surgeries, pregnant women, and patients with heart disease are not recommended to use.
7. If you feel any pain or abnormal skin reaction, please stop using immediately.

8. The cleaning and maintenance should be done when the power is off.
9. Do not store this product in a bathroom or any other wet room.
10. Do not allow children to use this product.
11. Do not use this product near flammable items (alcohol, gasoline, nail polish remover, spray, etc.).
12. Do not use the product while it is charging.
13. Do not disassemble or repair the product by yourself. For maintenance, please contact the store where you purchased the product.
14. Please use the matching USB line to charge to ensure that the charging voltage meets the product's requirements.
15. Do not charge when it is in use.
16. Do not use when taking a shower.

Usage Prohibitions

The following people shall not use:

1. Allergic people, those with skin diseases, dermatitis, skin allergies, or severe sun damage.
2. People with skin damage, eczema, swelling, or other similar conditions.
3. Individuals with physical, sensory, or mental issues, or those lacking experience.

Product Introduction

The Blackhead Remover Instrument is a beauty device that features blue light, exfoliation, pore cleansing, blackhead suction, grease suction, and other functions. It effectively improves the skin, making pores smooth and delicate while giving the skin a healthy, radiant appearance. The device includes 5 different suction heads, each offering a unique function.

Head Replaceable Design

The head, which has different functions, can be changed according to the varying needs of cosmetic effects.

Adjustable Work Dimension

There are high, medium, and low grades, with the option to switch between them based on your needs.

Easily to Be Washed

Five types of beauty heads can be easily removed for cleaning. It includes a filter sponge that effectively prevents dead skin from entering the machine. The sponge can be repeatedly cleaned for repeated use.

Compact and Easy to Carry

The rechargeable design is compact, lightweight, and easy to carry.

Large Capacity and Easy to Recharge

It features a built-in rechargeable lithium battery of 500mAh, which can be charged via the portable charging connector (USB). It can be recharged through a computer, mobile device, or any device with a USB interface, allowing use both at home or outdoors.

Blue-Ray Nursing

The 415nm blue light has anti-inflammatory effects, which help control and treat acne.

Component Description



Beauty Suction Head

1. Big Round Hole Beauty Head (Strong Suction): This beauty head provides strong suction and is primarily used for removing blackheads,

grease, dust, and cosmetic residue. Due to its strong suction, it is recommended for individuals with normal to resilient skin. Use with caution on fragile and thin skin. Limit usage to once a week, for no more than 5 minutes at a time, and adjust according to your skin condition.

Note: It is not recommended for use around the eyes but can be used on rough, wrinkled skin on the face, body, or elbows.

2. Silicone Head (Low Suction - Can Massage the Skin): The silicone head is designed for sensitive skin or areas that require special care. It provides a gentler cleansing effect, reducing friction and minimizing skin damage.

Note: It is not recommended for use around the eyes but can be used on rough, wrinkled skin on the face, body, or elbows.

3. Oval Hole Beauty Head (With Low Suction): The oval hole design is ideal for daily skincare, especially for targeting wrinkles around the corners of the mouth and similar areas. For best results, apply suction, quickly pull it away, and then apply suction again. You can also use it in an upward motion to lift the skin.

Note: It is suitable for use on the cheeks, eyes, nose, mouth, and other parts of the face.

4. Blue Ray Massage Head: The blue-ray head has a bactericidal effect that helps sterilize and purify the skin. This unique design also provides a gentle massage to reduce inflammation, control blackheads, and deeply nourish the skin. It is ideal for use on the cheeks and nose.



1



2



3



4

Replacement of Cosmetic Suction Head and Sponge

Remove the current suction head or sponge.

Insert the new suction head or sponge securely into place.

Ensure it is properly attached before use.



Specifications

Model No	GL-TM51
Material	ABS + PC
Battery	500mAh
Power	2.2W
Vacuum Pressure	≤ 60 kPa
Suction Levels	3 Levels
Noise Level	≤ 75 dB
Charging Port	USB-C
Charging Time	90 min
Usage Time	1.5 Hrs
Product Weight	132 g
Product Dimensions	40 x 43 x 179 mm

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

This **Green Lion** product comes with a 12-month warranty when purchased from our website, official store, or approved sellers.

To extend the warranty, please visit <https://www.greenlion.net/warranty> and fill out the form with your product details. Make sure to upload a clear picture of the product.

After we review and approve your request, you will receive a confirmation email stating that your product's warranty has been extended.

For more info, please check:

<https://www.greenlion.net/warranty>

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

Website: <https://www.greenlion.net/>

Sales: sales@greenlion.net

Service Support: Support@greenlion.net

Telegram: [@greenlion_global](https://t.me/greenlion_global)

Instagram: [@greenlion_global](https://www.instagram.com/greenlion_global)