

PORODO
LIFESTYLE



Porodo Lifestyle

**FACE & NECK Light Therapy
Massager Skin Rejuvenator**

SKU: PDHY48WH

Table of Contents

Product Overview	2
Specifications	2
Product Features	3
Product Design	3
Maintenance Guide	4
Safety Precautions	4
FCC Statement	5
Disposal	6
Warranty	6
Contact Us	6

Product Overview

The Porodo Lifestyle Face & Neck Light Therapy Massager is a skin rejuvenator designed to enhance your skincare routine with multiple features. It offers Ion Deep Cleaning for a thorough cleanse, Skin Firming Vibration to tone and tighten the skin, and Warm Massage Therapy for a soothing experience. Additionally, the massager uses Light Therapy with three modes—Purple, Blue, and Red—each targeting different skin concerns, such as calming and soothing, improving skin tone, and reducing fine lines. With four adjustable speed levels, this easy-to-use device is the perfect tool for anyone looking to improve their skin health and appearance.


Specifications

Rated Power	5W
Input	DC 5V/1A
Battery Capacity	450mAh
Charging Interface	Type-C
Modes	3 Light Modes (Blue, Red, Purple)
Recharge Time	3H
Product Size	73 x 92 x 100mm
Product Weight	200g
Model Number	PDHY48

Product Features

- 1. Ion Deep Cleaning: Cleans deeply to remove impurities and refresh the skin.
- 2. Skin Firming Vibration: Tightens and firms the skin for a smoother, youthful appearance.
- 3. Warm Massage Therapy: Offers a soothing and relaxing massage to enhance skin vitality.
- 4. Easy to Use: Simple operation for convenient skincare.
- 5. Purple Light Mode: Calms and soothes the skin, improving circulation and reducing puffiness for a balanced complexion.
- 6. Blue Light Mode: Helps balance oily skin, tightens pores, and boosts collagen production for enhanced skin vitality.
- 7. Red Light Mode: Accelerates blood circulation, boosts collagen, and firms skin while reducing wrinkles and fine lines.
- 8. Rechargeable: 3-hour recharge time for long-lasting use.
- 9. Type-C Charging: Fast and efficient charging with a modern Type-C interface.

Product Design



1. Power on/off button:

Long Press to switch on/off

Short Press to switch mode

2. Gear key

3. Color light

4. DC charging port

5. Gear display

6. Battery display

7. Work indicator light

8. Mode 1

9. Mode 2

10. Mode 3

11. Mode 4

Maintenance Guide

1. Store the product in a dry, room-temperature environment or in a storage box.
2. Avoid overcharging to protect the battery life.
3. Clean the product after charging and store it in the box if not in use for an extended period.
4. Use alcohol or a diluted detergent to wipe the exterior of the product.
5. The device may become slightly warm during use, but this is normal.
6. Do not rinse the device directly with water.
7. If not used for a long time, charge the battery at least once every 3 months.

Safety Precautions

1. Do not use this product with implanted medical devices, such as pacemakers, as it may cause equipment malfunction or pose a risk.
2. Avoid using the product if your palms have skin conditions like wounds, eczema, blisters, or peeling.
3. Do not use the product in high-humidity areas, such as bathrooms.
4. Use the product only for its intended purpose.
5. Do not disassemble, repair, or modify the product without authorization, as it may cause malfunction or injury.

6. Keep this product out of reach of children under six years old to prevent personal injury.

FCC Statement

1. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference; and

(2) This device must accept any interference received, including interference that may cause undesired operation.

2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for assistance.

RF warning statement:

The device has been evaluated to meet general RF exposure requirements. The device can be used in portable exposure conditions without restriction.

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

Products that you buy directly from our **Porodo** website or shop come with a 24-month warranty.

When you buy **Porodo** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at **porodo.net/warranty** and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:
porodo.net/warranty

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

info@porodo.net

Website: **porodo.net**

Service Support: **support@porodo.net**

Instagram: **[porodo](https://www.instagram.com/porodo)**