



Porodo Lifestyle

Smart Sonic Electric Toothbrush

SKU:PDLFSTETBWH

Table of Contents

Product Overview	2
Specifications	2
Intended Purpose	2
Care and Storage	3
Operations Instructions	3
Package Contents	4
FCC Statement	5
Disposal	6
Warranty	6
Contact Us	6

Product Overview

The Porodo Lifestyle Smart Sonic Electric Toothbrush offers advanced oral care with five customizable cleaning modes: Clean, Polish, White, Milder, and Massager. Its IPX7 water-resistant design ensures durability, while the smart timer helps maintain the ideal brushing time. With up to 40 days of use on a full charge, this toothbrush combines convenience and reliability. The ergonomic handle provides a comfortable grip, and the magnetic attachment allows for easy brush head replacement, making it a practical solution for effective oral hygiene.

Specifications

Input	5V/1A
Battery Capacity	300mAh
Recharge Time	2H
Working Time	Up to 40 Days
Water Resistance Level	IPX7
Material	ABS Body, Dupont Bristles
Product Size	19 × 23 ×245mm
Model Number	PDLFSTETBWH

Intended Purpose

This product is designed solely for personal oral hygiene use. It is not intended for commercial use, such as by medical or veterinary establishments. To ensure safe operation, please follow all instructions provided.

Care and Storage

- Clean the device according to the provided instructions.
- 2. Store the device in a dry, well-ventilated area, away from children.
- Protect the device from vibration and shocks during handling.
- **4.**The magnetic attachment works only at the designated spot on the toothbrush.



Operations Instructions



- 1. Charging
- a. Make sure the device is off.



- b. Connect the Type-C cable to the device.
- 0
- c. The LED indicator will flash slowly while charging.
- **d.** Unplug the cable when the battery is fully charged, indicated by a steady LED light.
- 2. Installing / Replacing the Head
- a. Ensure the device is off.
- b. Remove the head by pulling it off.
- c. Attach a new head by pushing it on and tightening it.
- 3. Program Selection
- a. Clean Standard cleaning mode (top LED light).
- b. Polish For teeth brightening and polishing (second LED light).
- c. White Removes discoloration, brightens, and polishes teeth (third LED light).
- d. Milder Gentle mode for sensitive teeth and gums (fourth LED light).
- e. Massager For gum massaging after brushing (bottom LED light).

- 4. Operation Mode and Timer
- **a.** To change the operation mode, press the on/off/toggle switch quickly until the desired program is selected.
- **b.** The brush features a built-in timer. After every 30 seconds, it pauses briefly to signal the elapsed time. A standard cleaning session lasts for 2 minutes, consisting of 4 intervals.
- **c.** Hold the button for 2 seconds to turn off. Your last mode is saved and will resume at the next use.
- 5. Brushing Instructions
- a. Wet the toothbrush head.
- b. Apply toothpaste to the brush.
- **c.** Insert the brush into your mouth before turning the device on to avoid splashing.
- d. Turn on the device by pressing the On/Off button.
- e. Start brushing from the gums, moving from left to right.
- f. Brush all teeth, starting from the gum area.
- g. Brush the inner sides of your teeth and repeat.
- h. Scrub the chewing surfaces carefully.
- i. After brushing, turn off the device and rinse the head under running water.

Notes:At first, your gums may bleed slightly due to the new type of irritation. This is normal, but if the condition persists for more than one week, consult your dentist.

Package Contents

Qty	Item
1	Brush Handle
2	Brush Heads
1	Magnetic Holder
1	USB Cable
1	User's Manual

Open the box and check that all items are included and in good condition. If any parts are damaged or missing, contact the dealer and do not use the product. Keep or recycle the packaging according to local regulations.

Note: Keep packaging materials (plastic bags, boxes, styrofoam, etc.) out of reach of children to prevent suffocation.

FCC Statement

- 1. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
- (1) This device may not cause harmful interference; and
- (2) This device must accept any interference received, including interference that may cause undesired operation.
- 2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for assistance.

RF warning statement:

The device has been evaluated to meet general RF exposure requirements. The device can be used in portable exposure conditions without restriction.

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

Products that you buy directly from our **Porodo** website or shop come with a 24-month warranty.

When you buy **Porodo** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at **https://www.porodo.net/warranty** and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:

https://www.porodo.net/warranty

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

info@porodo.net

Website: https://www.porodo.net/ Service Support: support@porodo.net

Instagram: porodo