

PORODO
LIFESTYLE



Porodo Lifestyle
AURA GLOW Light Therapy Mask
Wireless with 7-Color Cold Light

SKU: PDLFST206BK

Table of Contents

Specifications	2
Product Overview	2
Package Contents	3
Light Sources and functions	3
Frequently Asked Questions	5
Safety Precautions	6
Precautions Before Daily Use	8
Precautions for Daily Use	8
Precautions After Daily Use	9
Different Conditions and Targeted Solutions	9
Exclusions	10
Troubleshooting	11
Usage Instructions	12
Product Design	13
Disposal	14
Warranty	14
Contact Us	14

Specifications

Model	LFS206
Material	ABS + LEAD Beads
Input	5V/1A
Battery Capacity	400mAh
Charging Interface	Type-C
Working Time	Up to 50 Minutes
Rated Power	5W
Net Weight	298g
Product Size	83 x 192 x 230mm

Product Overview

This home photon beauty device utilizes multiple wavelengths of spectral light, designed to address a range of skin concerns at varying depths. It is widely recognized and trusted by users for its effectiveness and safety. The device works by emitting bio-active cold light, which does not generate high levels of heat. Its technology harnesses low-energy, pure light and unique modulation

techniques, making it gentle and suitable for sensitive skin. This innovative approach ensures minimal risk, with no known side effects, offering a natural, non-invasive solution for skin rejuvenation and improvement. For your safety, it is important to carefully read this instruction manual before using the product. After reading, please keep the manual in a safe place for future reference. If you decide to transfer the product to someone else, ensure that you also provide them with this manual.

Package Contents

Light therapy repair mask

Main body

Box

Connection cable

Light Sources and Functions

This product operates on the principle of light therapy, utilizing seven different colors of cold light, including red, blue, green, yellow, purple, cyan, and laser light. Each color offers distinct benefits, targeting various skin concerns and enhancing overall skin health. This advanced technology helps improve skin appearance by stimulating cellular activity, promoting collagen production, and providing rejuvenating effects.



Red Light (630nm)

Used for its soothing properties, red light therapy can promote circulation and help with skin rejuvenation.



Blue Light (460nm)

Known for its antibacterial effects, blue light can assist in the treatment of acne and skin irritation.



Green Light (515nm)

Green light is often used for reducing pigmentation and promoting skin healing, leaving the skin feeling calm and refreshed.



Violet Light (460nm)

Violet light offers anti-inflammatory benefits and is useful for skin detoxification and improving overall complexion.



Yellow Light (630nm)

Yellow light therapy can enhance skin elasticity, making it beneficial for anti-aging purposes and promoting a youthful appearance.



Cyan Light (460nm)

Cyan light is often employed for its cooling and calming effects, helping to balance the skin tone and reduce redness.



Laser Light (460nm)

Laser light provides more targeted therapy for deeper skin layers, aiding in the reduction of blemishes and other skin imperfections.

Frequently Asked Questions

Q1. How long does it take for a full treatment with the phototherapy mask?

The treatment duration for each mode is 15 minutes. After the treatment ends, the machine automatically shuts down. If you switch between modes, the machine will reset the timer. To ensure optimal use, the total cumulative exposure should not exceed 30 minutes per day.

Q2. Can I use the light therapy mask every day?

Yes, the light therapy mask can be used daily. However, it is recommended to choose the mode that best suits your skin condition on different days to achieve the most effective results.

Q3. Can I use the light therapy mask with a whitening essence?

Yes, you can. For the best results, it is advised to use the whitening essence after the light therapy treatment. Applying it before the session may reduce its effectiveness, so it's best to apply it afterward.

Q4. Can I use the phototherapy Mask with makeup on?

To achieve the best results, it is recommended to cleanse your face thoroughly and remove any makeup before using the phototherapy mask. This allows the light to penetrate effectively and enhances the treatment experience.

Q5. Can I use the mask while applying another mask?

If you want to use the phototherapy mask along with another face mask, we recommend using a medical-grade mask that is thin and contains simple ingredients. If the other mask is regular, it's better to apply it after the light therapy treatment to avoid interference with the results.

Q6. Who should avoid using the phototherapy mask?

The phototherapy mask is not recommended for children, pregnant women, or breastfeeding mothers. Individuals with photosensitive skin or specific skin conditions should consult with a doctor before using the mask to ensure it's safe for them.

Q7. Will I experience any skin reactions after using the phototherapy mask?

It is normal for the skin to feel slightly warm or dry after using LED light therapy. This occurs as the light stimulates skin cells. To alleviate this, apply a moisturizing mask or skin care products to hydrate your skin post-treatment.

Safety Precautions

- 1.** The product should be used in an environment where the temperature is between 10°C and 30°C to ensure optimal performance.
- 2.** Ensure that all buttons and indicators are functioning correctly before using the product for the first time or after it has been unused for a long period.
- 3.** If you have sensitive skin, particularly to light, it is advisable to perform a light sensitivity test before first use. If you experience any discomfort, stop using the product immediately.
- 4.** This product contains small parts. Keep it out of the reach of children to avoid accidental ingestion.
- 5.** Do not use this product for non-grooming purposes, as it may cause injury or malfunction.
- 6.** Do not use this product as a pet care tool, as it may cause harm to animals.

- 7.** This product should not be used by physically or mentally disabled persons without adult supervision, as improper use may lead to accidents or injury.
- 8.** Do not operate the product if it is malfunctioning. Using it in such a state may cause personal injury and damage to the device.
- 9.** Avoid using flammable substances (e.g., alcohol, gasoline, nail polish remover, or styling products) near the device. These substances may lead to electric shock, fire, or bodily injury.
- 10.** Never charge, use, or leave the device in a fire, heater, or hot location, as this may cause overheating, fire, or explosion. Always charge the device in a safe environment.
- 11.** Do not use the device near water, such as in bathrooms or during bathing, to prevent electric shock, short circuits, or other hazards.
- 12.** Avoid dropping or bumping the device. Mishandling the device may lead to malfunction due to internal damage from collisions.
- 13.** Do not use solvents (e.g., alcohol, gasoline, or nail polish remover) to clean the device, as they can damage its surface and functionality.
- 14.** Turn off the device when not in use, especially during cleaning, to avoid electric shock, short circuits, or fire risks.
- 15.** Ensure the power adapter used is from the manufacturer and matches the device's voltage requirements. Do not use extension cords, as they can cause electric shock or malfunction.
- 16.** Do not attempt to repair or modify the device yourself. Doing so may lead to electric shock, injury, or further damage. Contact the after-sales service for repair.
- 17.** If unusual smoke, noise, or odor occurs during operation, immediately disconnect the power adapter and contact after-sales service.
- 18.** Do not use the device with cosmetics or sprays that contain flammable ingredients, as this could cause electric shock, fire, or bodily injury.
- 19.** Do not use the device near electronic medical equipment such as pacemakers or ECG scanners, as it could interfere with medical devices and cause malfunctions.

Warnings:

- a. Failing to follow the directions for use may result in adverse effects or potential harm.
- b. Prolonged or repeated use of the product on the same area of skin may cause irritation or damage.
- c. Do not bypass or disable any safety features built into the product.

Precautions Before Daily Use

1. Individuals with allergies, skin conditions (e.g., eczema, dermatitis), or particularly sensitive skin may experience discomfort or adverse effects when using this product. Should you experience any negative effects, discontinue use immediately.
2. Avoid using the product on areas of the skin that are compromised (e.g., burns, cuts, blisters, or active infections). Do not use on areas where makeup or skin care products are still present. Clean the skin thoroughly before applying the product for skincare purposes.
3. Always remove any metal objects (e.g., jewelry, eyeglasses, or contact lenses) from the area to be treated before use.
4. The product (and adapter) should not be submerged in water under any circumstances.
5. Do not store the product in humid environments to prevent damage. Ensure the product is kept in a dry and cool place when not in use.
6. This product is not designed for use on children. Please keep out of reach of minors.

Precautions in Daily Use

1. Avoid direct eye exposure to the LED light for prolonged periods, as it may cause eye damage. If you notice any skin abnormalities, immediately stop using the product and seek assistance from after-sales support or consult a healthcare professional.

2. Select the appropriate mode based on your individual needs when using this product to ensure it suits your personal situation.
3. Limit exposure to the same area of skin to no more than 30 minutes per day to prevent potential skin damage.

Precautions After Daily Use

1. Ensure that both the main unit and adapter are kept dry and clean after each use.
2. When cleaning the product, avoid getting the connection line wet. After cleaning, make sure the product is fully dry before using it again.
3. A slight warmth and dryness on the skin after use is normal. For added comfort, use skincare products with moisturizing properties.
4. If you experience any discomfort or abnormal sensations while using the product, immediately stop use and consult a healthcare professional.

Different Conditions and Targeted Solutions

1. Managing Oily Skin & Preventing Acne Breakouts

Start by moisturizing with water and a mild milk-based product to keep skin hydrated without overloading it. Switch to "Acne Control Mode" to help reduce excess oil and prevent acne flare-ups. Consistency is key—ensure you follow a good skincare routine throughout this period for the best results.

2. Relieving Seasonal Dryness & Allergies (Beauty Recovery Period)

Apply gentle, non-irritating skincare products that are mild on the skin. Adjust to "Soothing Mode" to calm sensitive or reactive skin. Don't forget to protect your skin with sunscreen daily for the next 7 days to avoid further irritation from sun exposure.

3. Brightening and Rejuvenating Late-Night Skin

Start with a gentle cleansing routine to remove impurities. Switch to

"Brightening Mode" to revitalize tired skin and even out tone. To maximize the brightening effects, follow up with a whitening essence after treatment to nourish and restore your skin's natural glow.

Exclusions

This product should not be used by individuals with any of the following conditions, as it may lead to accidents or cause adverse skin and body reactions:

- 1. Light Sensitivity:** Illumination is not recommended for individuals with extensive skin diseases or those who are sensitive to light, including those with light allergies.
- 2. Medical Conditions:** Individuals with severe heart, liver, or blood system disorders, nervous system diseases, malignant tumors, facial neuralgia, hemophilia, febrile diseases, or conditions such as epilepsy that are triggered by light exposure should refrain from using this product.
- 3. Pregnancy:** Pregnant women are advised not to use this product, especially on the abdominal area. If necessary, consult a healthcare professional before use.
- 4. Medication Use:** If you are taking medication that may cause photosensitivity, please use this device with caution and strictly follow your healthcare provider's guidance.
- 5. Irritation or Discomfort:** If, during or after use, the skin becomes irritated (such as itching, swelling, redness, rashes, or blisters), or if discomfort is caused by light exposure, immediately stop using the product and consult a healthcare professional.

If at any point you feel discomfort during the irradiation process, stop the operation immediately and seek medical advice.

Troubleshooting

Concern	Presumed Cause	Prescription
Difficulty in powering on the device	The battery may be dead.	Recharge the battery and try turning it on.
Automatic shutdown of the equipment	+The device has reached the end of its service life.	+ After a 15-minute treatment, there might be a pause in the operation.
	+The device has been suspended for more than 10 minutes.	+ Pause for more than 10 minutes to allow for automatic shutdown.
Sensation of warmth on the face	Exposure to phototherapy could be a factor.	It is normal.
Dry skin sensation	Skin dehydration following phototherapy.	Use moisturizing skin care products.

Usage Instructions

- 1. Cleanse Your Face:** Begin by cleansing your face and allowing the water to air dry naturally. Alternatively, use a soft, absorbent towel to gently pat your face dry. Avoid rubbing your face with a coarse towel, as this can lead to enlarged pores and may damage dehydrated skin.
- 2. Prepare the Light Therapy Mask:** Take the light therapy repair mask and open both sides of the mask bracket. Position the mask under your forehead, ensuring it fits comfortably. Press the power button on the mask for about 1.5 seconds until it activates. You'll hear a sound indicating the switch to the red-light mode. After that, gently tap the light therapy mask to cycle through different light modes for skin repair. Be cautious not to look directly at the light source for prolonged periods during use.
- 3. Using the Photon Mask with Skincare Products:** After cleansing your face, apply your preferred skincare products or a disposable mask evenly onto your skin. Place the photon mask over your face, ensuring it stays in place while the red-light therapy is active. Let the mask remain on your skin for the recommended duration (typically 30-15 minutes). The power supply control box will automatically turn off the LED photon after the set time. Once the mask is removed, it's important to clean your face thoroughly. This will help maintain the effectiveness for future use.
- 4. Enhanced Results with a Translucent Silk Mask:** For even better results, consider using a translucent silk mask. This enhances the absorption of nutrients from your skincare products, especially when combined with the light therapy. The mask helps to lock in moisture, maximizing the benefits of the treatment.



1. Type-C Charging Port

2. +/- Keys and Universal Power Button

1. Charging: On the underside of the mask, you'll find the Type-C interface for charging.

Note: It will take approximately 2 hours to fully charge the mask.

2. Buttons Overview:

2.1 The + and - keys are used for adjusting settings such as intensity or mode.

2.2 The universal power button is used to power the mask on or off.

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

Products that you buy directly from our **Porodo** website or shop come with a 24-month warranty.

When you buy **Porodo** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at **porodo.net/warranty** and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:
porodo.net/warranty

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

info@porodo.net

Website: **porodo.net**

Service Support: **support@porodo.net**

Instagram: **porodo**